

CORMED

medizintechnik

Exclusive Distribution Partner





PRODUCT CATALOGUE

Extracorporeal Circuit Cannulas

SleekFlow® CV-S

Venous Return Cannula - Straight Plastic Tip





Ordering information

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
VC-PSR-12	12 Fr	4.0	380	1/4
VC-PSR-14	14 Fr	4.7	380	1/4
VC-PSR-16	16 Fr	5.3	380	1/4
VC-PSR-18	18 Fr	6.0	400	5/16
VC-PSR-20	20 Fr	6.7	400	5/16
VC-PSR-22	22 Fr	7.3	400	5/16
VC-PSR-24	24 Fr	8.0	400	5/16
VC-PSR-26	26 Fr	8.7	400	5/16
VC-PSR-28	28 Fr	9.3	400	3/8
VC-PSR-30	30 Fr	10.0	400	3/8
VC-PSR-32	32 Fr	10.7	400	3/8
VC-PSR-34	34 Fr	11.3	400	3/8
VC-PSR-36	36 Fr	12.0	400	1/2
VC-PSR-38	38 Fr	12.7	400	1/2
VC-PSR-40	40 Fr	13.3	400	1/2
VC-PSR-42	42 Fr	14.0	400	1/2

Extracorporeal Circuit Cannulas

SleekFlow® CV-S

Venous Return Cannula - Angle Plastic Tip




Ordering information

Order No.	Size	OD (mm)	Length (Mm)	Connector Port (inch)
VC-PCR-12	12 Fr	4.0	350	1/4
VC-PCR-14	14 Fr	4.7	350	1/4
VC-PCR-16	16 Fr	5.3	350	1/4
VC-PCR-18	18 Fr	6.0	370	5/16
VC-PCR-20	20 Fr	6.7	370	5/16
VC-PCR-22	22 Fr	7.3	370	5/16
VC-PCR-24	24 Fr	8.0	370	5/16
VC-PCR-26	26 Fr	8.7	370	5/16
VC-PCR-28	28 Fr	9.3	355	3/8
VC-PCR-30	30 Fr	10.0	355	3/8
VC-PCR-32	32 Fr	10.7	355	3/8
VC-PCR-34	34 Fr	11.3	355	3/8
VC-PCR-36	36 Fr	12.0	355	1/2
VC-PCR-38	38 Fr	12.7	355	1/2
VC-PCR-40	40 Fr	13.3	355	1/2
VC-PCR-42	42 Fr	14.0	355	1/2

SleekFlow® CV-S

Venous Return Cannula - Angle Metal Tip





Ordering information

Order No.	Size	Od (mm)	Length (mm)	Connector Port (inch)
VC-MCR-12	12 Fr	4.0	395	1/4
VC-MCR-14	14 Fr	4.7	395	1/4
VC-MCR-16	16 Fr	5.3	395	1/4
VC-MCR-18	18 Fr	6.0	425	5/16
VC-MCR-20	20 Fr	6.7	425	5/16
VC-MCR-22	22 Fr	7.3	425	5/16
VC-MCR-24	24 Fr	8.0	425	5/16
VC-MCR-26	26 Fr	8.7	425	5/16
VC-MCR-28	28 Fr	9.3	425	3/8

Extracorporeal Circuit Cannulas

SleekFlow® CV-D

Venous Return Cannula - Two Stage Venous Cannula Round Lumen





Ordering information

Order No.	Size	1 st Stage (mm)	2 nd Stage (mm)	Length (mm)	Connector Port (inch)
		OD	OD		
VC-TSR-2428	24/28 Fr	8.0	9.3	390	3/8
VC-TSR-2836	28/36 Fr	9.3	12.0	390	1/2
VC-TSR-3038	30/38 Fr	10.0	12.7	390	1/2
VC-TSR-3238	32/38 Fr	10.7	12.7	390	1/2
VC-TSR-3438	34/38 Fr	11.3	12.7	390	1/2
VC-TSR-3440	34/40 Fr	11.3	13.3	390	1/2
VC-TSR-3446	34/46 Fr	11.3	15.3	400	1/2
VC-TSR-3448	34/48 Fr	11.3	16.0	400	1/2
VC-TSR-3646	36/46 Fr	12.0	15.3	400	1/2
VC-TSR-3650	36/50 Fr	12.0	16.0	400	1/2
VC-TSR-3850	38/50 Fr	12.7	16.0	400	1/2

Extracorporeal Circuit Cannulas

SleekFlow® CV-D

Venous Return Cannula - Two Stage Venous Cannula Oval Lumen




Ordering information

Order No.	Size	1 st Stage (mm)	2 nd Stage (mm)	Length (mm)	Connector Port (inch)
		OD	OD		
VC-TSO-2428	24/28 Fr	8.0	9.3	390	3/8
VC-TSO-2836	28/36 Fr	9.3	12.0	390	1/2
VC-TSO-3038	30/38 Fr	10.0	12.7	390	1/2
VC-TSO-3238	32/38 Fr	10.7	12.7	390	1/2
VC-TSO-3438	34/38 Fr	11.3	12.7	390	1/2
VC-TSO-3440	34/40 Fr	11.3	13.3	390	1/2
VC-TSO-3446	34/46 Fr	11.3	15.3	400	1/2
VC-TSO-3448	34/48 Fr	11.3	16.0	400	1/2
VC-TSO-3646	36/46 Fr	12.0	15.3	400	1/2
VC-TSO-3650	36/50 Fr	12.0	16.0	400	1/2
VC-TSO-3850	38/50 Fr	12.7	16.0	400	1/2

SleekFlow[®] AO-Plus

Arterial Cannula Straight Plastic Tip, Wirewound Body, with Introducer





Ordering information

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PSR-ST-08	8 Fr	2.7	260	1/4
IFAC-PSR-ST-10	10 Fr	3.3	260	1/4
IFAC-PSR-ST-12	12 Fr	4.0	260	1/4
IFAC-PSR-ST-14	14 Fr	4.7	260	1/4
IFAC-PSR-ST-16	16 Fr	5.3	260	1/4

SleekFlow[®] AO

Arterial Cannula Straight Plastic Tip, Wirewound Body




Ordering information (Without Suture Ring)

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PSR-08	8 Fr	2.7	211	1/4
IFAC-PSR-10	10 Fr	3.3	211	1/4
IFAC-PSR-12	12 Fr	4.0	211	1/4
IFAC-PSR-14	14 Fr	4.7	211	1/4
IFAC-PSR-16	16 Fr	5.3	211	1/4
IFAC-PSR-18	18 Fr	6.0	255	3/8
IFAC-PSR-20	20 Fr	6.7	255	3/8
IFAC-PSR-22	22 Fr	7.3	255	3/8
IFAC-PSR-24	24 Fr	8.0	255	3/8

SleekFlow[®] AO

Arterial Cannula Straight Plastic Tip, Wirewound Body, with Suture Ring





Ordering information (With Suture Ring)

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PSR-SR-14	14 Fr	4.7	211	1/4
IFAC-PSR-SR-16	16 Fr	5.3	211	1/4
IFAC-PSR-SR-18	18 Fr	6.0	255	3/8
IFAC-PSR-SR-20	20 Fr	6.7	255	3/8
IFAC-PSR-SR-22	22 Fr	7.3	255	3/8
IFAC-PSR-SR-24	24 Fr	8.0	255	3/8

SleekFlow® AO

Arterial Cannula Straight Plastic Tip, not Wirewound Body





Ordering information (Without Suture Ring)

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PSN-08	8 Fr	2.7	211	1/4
IFAC-PSN-10	10 Fr	3.3	211	1/4
IFAC-PSN-12	12 Fr	4.0	211	1/4
IFAC-PSN-14	14 Fr	4.7	211	1/4
IFAC-PSN-16	16 Fr	5.3	211	1/4
IFAC-PSN-18	18 Fr	6.0	255	3/8
IFAC-PSN-20	20 Fr	6.7	255	3/8
IFAC-PSN-22	22 Fr	7.3	255	3/8
IFAC-PSN-24	24 Fr	8.0	255	3/8

Extracorporeal Circuit Cannulas

SleekFlow[®] AO

Arterial Cannula Straight Plastic Tip, not Wirewound Body, with Suture Ring





Ordering information (With Suture Ring)

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PSN-SR-14	14 Fr	4.7	211	1/4
IFAC-PSN-SR-16	16 Fr	5.3	211	1/4
IFAC-PSN-SR-18	18 Fr	6.0	255	3/8
IFAC-PSN-SR-20	20 Fr	6.7	255	3/8
IFAC-PSN-SR-22	22 Fr	7.3	255	3/8
IFAC-PSN-SR-24	24 Fr	8.0	255	3/8

SleekFlow[®] AO

Arterial Cannula Angle Plastic Tip, Wirewound Body




Ordering information

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PCR-12	12 Fr	4.0	223	1/4
IFAC-PCR-14	14 Fr	4.7	223	1/4
IFAC-PCR-16	16 Fr	5.3	223	1/4
IFAC-PCR-18	18 Fr	6.0	264	3/8
IFAC-PCR-20	20 Fr	6.7	264	3/8
IFAC-PCR-22	22 Fr	7.3	264	3/8
IFAC-PCR-24	24 Fr	8.0	264	3/8

SleekFlow[®] AO

Arterial Cannula Angle Plastic Tip, not Wirewound Body



Ordering information

Order No.	Size	OD (mm)	Length (mm)	Connector Port (inch)
IFAC-PCR-12	12 Fr	4.0	223	1/4
IFAC-PCR-14	14 Fr	4.7	223	1/4
IFAC-PCR-16	16 Fr	5.3	223	1/4
IFAC-PCR-18	18 Fr	6.0	264	3/8
IFAC-PCR-20	20 Fr	6.7	264	3/8
IFAC-PCR-22	22 Fr	7.3	264	3/8
IFAC-PCR-24	24 Fr	8.0	264	3/8

Cardioplegia Circuit Cannulas

LifeWell® TUB


Cardioplegia Circuit Cannula - Aortic Root Cannula

Medical grade PVC material for the body.

The special stainless steel treatment makes the needle rigid and tenacious.

Molded tip resists kinking.


With Vent Line



Ordering information

Order No.	Description
ARC-Y-15	With vent line, 15Ga, fixation flange, female luer connector
ARC-Y-16	With vent line, 16Ga, fixation flange, female luer connector

Without Vent Line



Ordering information

Order No.	Description
ARC-S-15	Without vent line, 15Ga, fixation flange, female luer connector
ARC-S-16	Without vent line, 16Ga, fixation flange, female luer connector


LifeWell® TUB

Cardioplegia Circuit Cannula - Coronary Artery Ostial Cannula 45° / 90°

These cannulas feature flanged, radiopaque basket tips.

All cannulas have female luer connection site.


45°



Ordering information

Order No.	Description
CAOC-45-10	45 degree, 10Fr (3.3mm)
CAOC-45-12	45 degree, 12Fr (4.0mm)
CAOC-45-14	45 degree, 14Fr (4.7mm)

90°



Ordering information

Order No.	Description
CAOC-90-10	90 degree, 10Fr (3.3mm)
CAOC-90-12	90 degree, 12Fr (4.0mm)
CAOC-90-14	90 degree, 14Fr (4.7mm)



Extracorporeal Circulation Series

Suction Circuit Cannulas

Suction Circuit Cannulas are intended to provide suction within and around the heart.

Suction Circuit Cannulae

SleekFlow® CC-LT

Vent Catheter Type A / Type B

Type A

Vent Catheter is designed to drain blood from left atrium and left ventricle.



Ordering information

Order No.	Size	OD of Stylet (mm)	ID of Tube (mm)	OD of Tube (mm)	Length (mm)	Number Of Side Holes	Connector Pore (inch)	Remark
SCA-18	18 Fr	2.8	3.5	5.8	358	5	1/4	With stylet
SCA-14	14 Fr	2.1	2.5	4.5	316	5	1/4	With stylet
SCA-12	12 Fr	-	2.5	4.0	308	6	1/4	Without stylet

Type B

Malleable suction catheter. The tip of this cannula can be curved to any angle at doctor's wish.



Ordering information


Order No.	Size	ID of Tube (mm)	OD of Tube (mm)	Length (mm)	Number Of Side Holes	Connector Pore (inch)
SCB-18	18 Fr	2.5	5.8	350	5	1/4
SCB-14	14 Fr	2.5	4.5	310	6	1/4
SCB-12	12 Fr	2.5	4.0	305	6	1/4

SleekFlow® CC-RT

Rigid Suction Wand / Flexible Suction Wand

Rigid Suction Wand

Feature maximum flow and minimum trauma. A plastic grip is equipped for comfortable operation. The tip of the pediatric model has been made softer for gentle suction.



Ordering information

Order No.	Size	ID of Tube (mm)	OD of Tube (mm)	Length (mm)	Connector Pore (inch)
SCR-18	18 Fr	4.8	5.8	295	1/4
SCR-14	14 Fr	3.5	4.5	230	1/4

Flexible Suction Wand

Feature a soft tip to maximize the suction flow and minimize the tissue trauma. The tip design is ideal for atraumatic suction within the heart chambers.



Ordering information

Order No.	Size	ID of Tube (mm)	OD of Tube (mm)	Length (mm)	Connector Pore (inch)
SCF-18	18 Fr	4.0	6.0	440	1/4
SCF-14	14 Fr	3.0	5.0	380	1/4



Kontakt:

Technische Beratung & Vertrieb

André Ginther

Vertriebsleitung

0160/91730834

andre.ginther@cormed.de

Innendienst:

Anna Maria Schlensak

Legal and Compliance Analyst /

Qualitätsmanagement

02954 / 97 55 – 28

anna.schlensak@cormed.de

Elena Wiens

CustomerService

02954/9755 – 10

elena.wiens@cormed.de

Elena Gross

Customer Service

02954/9755 – 55

elena.gross@cormed.de